MEDIKIT™

MK-NCov V.21 / 2020

APRIL 2020

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COVID-19 Antibody Rapid Test Kit

Lateral-Flow Immunoassay Rapid Test

97.51%
Accuracy

For use by
medical
professionals
only

SCAN ME

The world's first and only COVID-19 Antibody test kit with App enabled tracking, result interpretation and data collection, aiding healthcare services fight the virus and flatten the curve.

AVAILABLE ON GOOGLE PLAY
AND THE APPSTORE

PRODUCT SPECIFICATIONS

Product Name:	COVID-19 Antibody Test
Certificated Available:	CE, MHRA, COA, ISO*, Clinical Assessment Data
Product Manufacture:	Quick Energy Technologies Ltd

MEDIKIT SINGLE TEST PACK:

**ONE YEAR
WARRANTY**

MEDIKIT PACK CONTENT:

- 1 x Test Cassette**
- 1 x Phosphate Buffer Saline**
- 1 x Pipette**
- 1 x Alcohol Wipe**
- 1 x Lancet**
- 1 x QR Code**

MEDIKIT MULTI-PACK:


MEDIKIT PACK CONTENT:

- 25 x Test Cassette**
- 1 x Phosphate Buffer Saline (4ml)**
- 1 x QR Code**

Download the Medikit App to help keep track of **COVID-19** cases to help national Health services fight the virus and flatten the curve.

- 1 Download the app and follow the steps on how to use the test.
- 2 Allow access to your location.
- 3 Enter your age and gender.
- 4 Prepare and complete the MediKit test kit.
- 5 Wait **10 minutes** – check results and use the apps scanner to photograph the cassette to help keep track of cases.

COMPETITOR ANALYSIS

	CE Certification	MHRA	European Product	Lancet Included	Single-test Packaging	24/7 Online Support	Mobile App
MEDIKIT	✓	✓	✓	✓	✓	✓	✓
 BIOZEK	✓	✓	✓	✗	✗	✗	✗
 Wondfo	✓	✓	✗	✗	✗	✗	✗
 innovita	✓	✓	✗	✗	✗	✗	✗

TEST DATA

The MediKit Covid-19 IgG/IgM Test was compared with a leading commercial PCR (Polymerase chain reaction).

MediKit Covid-19 IgG/IgM Test	Diagnostics Sample				Total	
	Positive		Negative			
Positive	A	129	B	0	A + B	129
Negative	C	13	D	381	C + D	391
Total	A + C	142	B + D	381	ABCD	523

Sensitivity = $A/(A+C)\% = 90.85\%$ (95% credibility interval: 84.85~95.04%)

Specificity = $D/(B+D)\% = 100.00\%$ (95% credibility interval: 99.04%~100%)






Total Accuracy = $(A+D)/(ABCD)\% = 97.51\%$ (95% credibility interval: 95.79%~98.67%)

	Values	Percentage (95% confidence interval)
positive coincidence rate	129/142	90.85% (84.8 ~ 95.04%)
negative coincidence rate	381/381	100.00% (99.04% ~ 100.00%)
positive predictive value	129/129	100.00% (97.18% ~ 100.00%)
negative predictive value	381/394	96.70% (94.42% ~ 98.23%)
total coincidence rate	510/523	97.51% (95.79% ~ 98.67%)

PACKING SHEET

Dimensions

CE CERTIFICATION

 **NEWSCEN**

CE DECLARATION OF CONFORMITY

According to the In Vitro Diagnostic Medical Devices
Directive 98/79/EC

.....

Manufacturer: **NewScen Coast Bio-Pharmaceutical Co., Ltd.**
Address: No. 65, 6th Street, Tianjin TEDA, 300457, Tianjin, PEOPLE'S REPUBLIC OF CHINA

European Representative: Wellkang Ltd
Address: 16 Castle St, Dover, Kent, CT16 1PW, England, UK

Product: COVID-19 IgG/IgM Rapid Test Cassette

Classification: Others

Conformity assessment route: Annex III, excluding 6

*We, the Manufacturer, herewith declare with sole responsibility that our product mentioned above meets the provisions of the **Directive 98/79/EC** of the European Parliament and of the Council on In Vitro Diagnostic Medical Devices.*

Applicable Standards: EN ISO 13485:2016 EN 13612:2002
 EN ISO 14971:2012 EN ISO 23640:2015
 EN ISO 18113-1:2011 EN ISO 18113-2:2011
 EN ISO 15223-1:2016 EN 13641:2002

Signed on 12th day Mar. month/ of 20 20
Place TEDA city Tianjin Country China

Signature (General Manager): [Signature]

Medicines & Healthcare products
Regulatory Agency

MHRA

MHRA

10 South Colonnade
Canary Wharf
London
E14 4PU
United Kingdom

Our Ref: IVD000904

Dr Edward Wang
Wellkang Ltd
Suite B
29 Harley Street
London
W1G 9QR

27 March 2020

Dear Dr Wang

IN VITRO DIAGNOSTIC MEDICAL DEVICES REGULATIONS 2002: REGULATION 44
Registration of manufacturers of In-Vitro Diagnostic Medical Devices
and devices for Performance Evaluation

Thank you for informing the Competent Authority of the change to the original notification dated (date the registration was registered); **Manufacturers Name:- Newscen Coast Bio-Pharmaceutical Co., Ltd.** located at **Manufacturers Address:- No.65, 6th Street, Tianjin TEDA, Tianjin, China 300457** for whom you are acting as the authorised representative and for supplying the medical device information.

The change(s) to your registration has been recorded based on your declaration that you have determined that the device(s) fall within the definition of "in vitro diagnostic medical device", and that you have classified it/them correctly taking into account the intended purpose(s) and mode(s) of action. In accepting your registration, I should make clear that the Competent Authority does not examine each individual notification and therefore cannot and does not necessarily endorse these determinations. Neither does this letter represent any form of accreditation or approval by the UK Competent Authority.

Your registration is based upon your declaration on the RG3 form and means that you should now be operating under the In Vitro Diagnostic Medical Devices Directive and the 2002 Regulations for the products you asked us to register, by fully complying with the essential requirements, CE marking those products or labelling them as such.

If you stop placing devices on the market or if you are not complying with the Regulations you should inform us as required by the Regulations. You should be aware that it is an offence to place on the market CE marked devices that do not comply with the Regulations.

The information you provided has been recorded against the reference number shown at the top of this letter, which we ask you to quote in all future correspondence and communications.

Please inform us of any changes to:

- the company information
- additional generic groups of devices or, for Annex II or Self-Test devices, additional devices
- discontinuation of a generic group of devices or, for Annex II or Self-Test devices, discontinuation of devices

Medicines & Healthcare products
Regulatory Agency

Please use RG3, the Registration form, to tell us about any of these changes. A fee of £70 is payable for each change or set of changes notified.

Thank you for registering the following generic groups of devices

1. **Part 5: IVDs which are not Annex II and not self-test devices**
- 2.
3. **For reagents, reagent products, calibration and control materials:**
4. **group by common technological characteristics and/or analytes**
- 5.
6. **New products:**
7. **None**
- 8.
9. **For performance evaluation:**
10. **None**
- 11.
12. **Neither:**
13. **IC&II Hardware + accessories + consumables + software**
14. **C-Reactive Protein**
15. **Procalcitonin**
16. **Alphafetoprotein**
17. **Carcinoembryonic Antigen**
18. **Cardiac Marker Controls**
19. **D-Dimer Rapid Test**
20. **Albumin (IC)** CC =>11.02.01.01
21. **Thyroid Stimulating Hormone**
22. **Coronavirus**
- 23.
- 24.
25. **For other IVDs, group by appropriate indications**
- 26.
27. **New products:**
28. **None**
- 29.
30. **For performance evaluation:**
31. **None**
- 32.
33. **Neither:**
34. **None**
- 35.
- 36.
37. **Part 6: IVDs which are Annex II or self-test devices**
- 38.
39. **For reagents, reagent products, calibration and control materials:**
40. **group by common technological characteristics and/or analytes**
- 41.
42. **New products:**
43. **None**
- 44.
45. **For performance evaluation:**
46. **None**
- 47.
48. **Neither:**
49. **None**
- 50.
- 51.
52. **For other IVDs, group by appropriate indications**

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Certificate

No. Q5 095791 0008 Rev. 01

Holder of Certificate: **NewScen Coast
Bio-Pharmaceutical Co., Ltd.**
No. 65, 6th Street
Tianjin TEDA
300457 Tianjin
PEOPLE'S REPUBLIC OF CHINA

Facility(ies): NewScen Coast Bio-Pharmaceutical Co., Ltd.
No. 65, 6th Street, Tianjin TEDA, 300457 Tianjin, PEOPLE'S
REPUBLIC OF CHINA

Certification Mark:

Scope of Certificate: **Design and Development,
Production and Distribution of
Lateral Flow Technique based
In Vitro Diagnostic Medical Devices
and Readers for Lateral-Flow Devices**

Applied Standard(s): EN ISO 13485:2016
Medical devices - Quality management systems -
Requirements for regulatory purposes
(ISO 13485:2016)
DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.: SH19107304

Valid from: 2019-10-21

Valid until: 2022-10-20

Date, 2019-09-19

Stefan Preiß
Head of Certification/Notified Body

America

CERTIFICATE

No. QS5 095791 0009 Rev. 01

Certificate Holder:

NewScen Coast
Bio-Pharmaceutical Co., Ltd.
No. 65, 6th Street
Tianjin TEDA
300457 Tianjin
PEOPLE'S REPUBLIC OF CHINA

Certification Mark:

Scope of Certificate:

Design and Development, Production and Distribution
of Lateral-Flow Technique based In-Vitro Diagnostic
Medical Devices and Readers for Lateral-Flow Devices

Standard(s):

ISO 9001:2015

The Certification Body of TÜV SÜD America Inc. certifies that the company mentioned above has established and is maintaining a quality management system that meets the requirements of the listed standards.

Report No.:

SH19107304

Effective Date:

2020-03-09

Expiry Date:

2022-11-24

Page 1 of 1

Date of Issue: 2020-03-10

(Tina Israel)
Manager, US Certification Body,
Medical and Health Services

The antibody-based tests: Report & analysis

Healthcare settings pose a significant threat to patients and staff due to the increased risk of disease transmission through airborne pathogens. In light of this, during the COVID19 pandemic, this has resulted in an all but complete cessation of elective and primary care provisions of medical and allied care professions.

By the very nature of routine, invasive procedures generating aerosols of high risk biological fluids, dental practice may offer (worst case scenario/best case to study)

This is, in part, contributed to by the intrinsic exposure to high risk biological fluids.

In routine practice this is mitigated through routine use of protective equipment, meticulous hygiene and disinfection measures. With the advent of the COVID19 pandemic, and in the absence of effective countermeasures and diagnostics, this has led to the complete cessation of routine primary and specialised dental care provisions across England.

Whilst this has, undoubtedly, reduced the number of nosocomial infections arising from these settings, this has resulted in significant morbidity arising from the absence of these care provisions. There is a substantial body of debate relating to a staged restarting of these healthcare provisions. One such approach proposed comprises a combination of measures to mitigate risks within practices and testing of patients & staff returning to work. Widespread implementation of antibody based tests may provide insights into current and historical exposures to SARS-nCoV2-19, reflective of these disease states and as such a combined, risk-based approach ultimately enabling a resumption of routine, dental care provisions.

Whilst an increasing, often conflicting, body of data is becoming available regarding point of care tests for COVID19, we aimed to explore the applicability, acceptance and potential implementation of such point-of-care tests within routine primary care practice, lacking specialist testing facilities available to hospitals and within a short period of time.

Secondary endpoints of the study were to assess the prevalence of current and historical exposure and factors influencing these within these specific settings.

We attempted, using a representative sample of such settings, to explore this with a view to testing 13 dental practices with NHS and private patient bases within the Greater London area were invited to this study. Within these practice clusters, consenting participants undertook voluntary and anonymised testing based on serum IgM and IgG antibodies against SARS-nCoV2-

19, based on a CE-marked, lateral-flow immunoprecipitation kit with colorimetric read-out (COVID19 Antibody Kit, MediKit, London, United Kingdom). In addition to anonymised test results, participants were asked to provide basic information on gender, age brackets, symptoms, exposure- and work-based risk classification, confirmatory testing and ease of use/acceptance. The tests were conducted between March 23rd (6,700 cases in the UK confirmed) and April 25th, 2020 (149,500 cases UK-wide). Pseudonymised cassettes were returned and findings confirmed by the investigators.

Returns and demographics

In total, 106 complete kits with instructions as well as questionnaires were distributed, with a median of 8 tests per practice cluster conducted (range 4-16), >95% of tests provided were conducted (n=101), the remainder returned (please check 100 tests were valid (99%), with one a failure likely to have been an inappropriate sample application, precluding either positive or negative read-out. 96 questionnaires were returned, giving a matching return rate of 95.8%.

The participants demographics are reported in the table below.


Factor	n=	Percentage	Comments
Gender			
Male	32	31.7	
Female	69	68.3	
Age range			
30-40	48	50.5	
40-50	38	40	
50-60	9	9.5	
<30 or >60	0	0	
High risk symptoms (1/2/ALL)	4/3/5	4.2/3.1/5.2	n=12 total (13.5%)

Temperature >37.8°C	11	11.5	
Shortness of breath	5	5.2	
Persistent Cough	9	9.4	
Low risk symptom (EITHER/ALL)	4/9	4.2/9.4	12/13 also had HR symptoms
Loss of taste/smell	12	12.5	1 case no HR symp
RT symptoms	10	10.4	
Aggravating MH	8	8.3	All had >1 both HR/LR symptom
Exposure risk factor (HIGH/LOW)	87/9	90.6/9.4	
HIGH RISK:HW/NHW	78/8	90.7	
Swab test (Done/+ve)	4/3	4.2/3.2	All >1 symptom HR/LR as well
Ease of use (Y/N)	96	100	
Ease of read out	95	99	Included initially failed as negative
Ease of interpretation	95	99	Included initially failed as negative

*Numeric discrepancies in total numbers due to incomplete forms/withheld information, percentages expressed as percentage of total for the variable examined.

General comments

Overall, the study demonstrated a significant bias in terms of sampling by gender towards females, likely explained by the predominance of female dentists, allied registrants and team members within a typical practice mix. Most participants defined themselves as within high risk groups (91%), with predominantly by healthcare workers contributing to this subgroup with 90.7% (n=78).

All participants stated that they were able carry out the test easily (one duplicate), found the instructions easy to understand (100%) and felt they were able to understand the differential meanings of test results (100%, n=101). This is reflected in an overall successful completion rate of 99%. No failures of controls were observed. Again, 100% of participants stated that they would consider paying for these tests in future for practice use.

Correlations across the data


8% of the participants disclosed a medical history putting them at elevated risk, and remarkably, all displayed >2 symptoms with at least one high-specificity symptoms of SARS-nCOV2-19 infection. This contrasted with 5 reports of one or more symptoms of either symptoms group for the remaining 88 participants without any medical history putting the participant at high risk of adverse COVID19 outcomes. This means 100% versus 5.6% for MH and no-MH groups, respectively.

In other terms, a positive risk MH meant a 15.6 times higher likelihood of COVID symptoms being reported.

Within the reporting period, 13 participants displayed one or more symptom consistent with any course of SARS-nCOV2-19 infection, of which 12 had symptoms comprising isolation criteria. The

single most common presenting symptom was anosmia/agusia (loss of smell/taste, n=12), followed by fever (n=11). All patients that underwent swab testing disclosed a history of COVID19 symptoms, including specific isolation criteria in all cases.

Antibody based testing

18% of tested individuals were positive overall by the tests used.

17% displayed IgG only, suggestive of past infection and likely immunity to SARS-nCoV19-2.

1% had a combined IgM/IgG signal.

Healthcare workers were over-represented and could be aligned into practice clusters.

Considering isolation was in place at the point of testing, this would account for only one IgM positive result, as would be expected through social distancing and supportive of exposure prior to the government lock down.


Of note were interesting trends for clusters in practices tested, suggesting past exposure, possibly within the healthcare setting.

	YES/NO (n)	Rate YES	Comment
Successful test*	100/1	99%	One failure (I would attribute to execution rather than due to the test)
Positive band (any)	18/83	18%	
IgM only	0	0%	
IgM&IgG	1	1%	
IgG only	17	17%	

*a

successful test was defined as having a clearly distinguishable control band at the completion of the test. One sample was flooded with blood, likely due to puncture of the test field during the test, precluding a read out. This test was repeated, revealing a valid, negative read out.

62% of receiving clinicians reported one or more cases, v 38% no cases.

The average and median rate of positive Antibody within practices where at least one positive case was present were 37.5 and 41.2 % respectively.

Correlations across the tests

For the IgM/IgG antibody EIA testing, a complete correlation was noted of IgG counts and positive swabs, suggesting no false negatives. This was limited by a relatively small of swab tested participants reflective of current testing rates (n=3 confirmed, 4 tested). The one case that had a negative swab test result but multiple high-risk symptoms, tested positive by IgG test suggesting indeed a positive history and in keeping with the reported false negative rate of approximately 25 % for swab tests with qRT-PCR testing methodologies.

Correlation of test positivity			
Agreement with Swab test indication**	4/4	100	Meaning all pt having had swab test did test +ve EIA Ab test
Agreement with positive swab test	3/3	100	See comment in text

Antibody positive cases

In addition to the cross-validation through the few available swab based nucleic acid tests, the specificity and sensitivity was highly correlated to the presence of COVID19 symptoms. No participant showing any 5 symptoms suggestive of COVID19 had a positive antibody test (n=83).

Of significant interest was the observation of 6 of 18 participants reporting no symptoms, suggesting that these were asymptomatic carriers, however these could predominantly be placed within clusters of practices. This may implicate a significant risk posed by untested staff.

All of the subjects that showed COVID 19 or respiratory/sensory symptoms were within the positive group, thereby supporting the validity of tests and underpinning the hypothesis of a large, undetected number of individuals exposed to COVID19, particularly within health settings.

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